|   | FR                                 | DE                       | NL                         | EN                         |
|---|------------------------------------|--------------------------|----------------------------|----------------------------|
| ① | Bouton Marche/Arrêt                | Ein-/Ausschalter         | Knop Start/Stop            | On/Off button              |
| ② | Bouton de selection du mode        | Modus-Taste              | Selectieknop werkwijze     | Mode selection button      |
| ③ | Bouton de sélection de l'intensité | Intensitäts-Taste        | Selectieknop intensiteit   | Intensity selection button |
| ④ | Bouton de réglage de la chaleur    | Wärmetaste               | Regelknop warmte           | Heat setting button        |
| ⑤ | Témoin lumineux Marche/Arrêt       | Leuchtanzeige Ein/Aus    | Controlelampje Start/Stop  | On/Off indicator light     |
| ⑥ | Témoin de charge                   | Ladezeige                | Laadlampje                 | Charge indicator light     |
| ⑦ | Témoin lumineux du mode            | Leuchtanzeige Modus      | Controlelampje modus       | Mode indicator light       |
| ⑧ | Témoin lumineux intensité          | Leuchtanzeige Intensität | Controlelampje intensiteit | Intensity indicator light  |
| ⑨ | Témoin lumineux chaleur            | Leuchtanzeige Wärme      | Controlelampje warmte      | Heat indicator light       |

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**Warning:** Please read this instruction manual carefully before using the instrument and keep it for future use.

#### **A. SAFETY INSTRUCTIONS**

1. The device is intended exclusively for private and not for commercial use. Only use this device as describe in this instruction manual.
2. Do not place any part of this appliance in water or any other liquid.
3. Do not handle the unit or plug with wet hands.
4. When disconnecting the cable, always pull at the plug; do not pull on the cable itself.
5. When using the device, always pull the power cord completely from its housing to prevent overheating.
6. If the cord is damaged, it must be replaced by the manufacturer or its service agent or a similarly qualified person in order to avoid an electric shock.
7. Never try to repair it by yourself. It must be repaired by the manufacturer or its service agent or a similarly qualified person in order to avoid an electric shock.
8. Do not use outdoors.
9. Do not let the power cord on the edge of a table or work surface, and prevent it comes into contact with hot surfaces.
10. Keep it away from heat sources such as radiators, to avoid deformation of plastic parts.
11. Do not use coarse or abrasive sponges/cloths or steel wool to clean the appliance.
12. This unit is not designed to be operated by means of an external timer or a separate remote-control system.
13. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capacities, or by persons lacking in experience or knowledge, unless they are supervised or instructed on the use of this apparatus by a person responsible for their safety. Children should be supervised to ensure they do not play with this appliance.

#### **B. SPECIFIC SAFETY INSTRUCTION**

1. Do not let water enter the product or be splashed near the surface of the product.
2. Do not operate the button with wet hands.
3. Keep the buttons from water or other liquid.
4. Water contact with electricity is dangerous. Do not leave the product in a humid environment (such as in the bathroom or shower room or swimming pool).
5. Before connecting the product to the power supply, check that the product voltage and supply voltage are the same.
6. Product is not intended for use by the physically or mentally impaired (including children), unless there is supervision or guidance by the people responsible for their safety when they use the product.
7. Children should not use this product. They should be under supervision to ensure that they will not treat this product as a toy.
8. If an unusual sound or noise is being heard from the product, turn off the product and immediately disconnect the power connection. Please contact Service center immediately.
9. Do not use this product in high temperature conditions or beside hot objects such as heaters, stoves.
10. Do not apply medical or cosmetic creams on the product or on your foot before using this product.
11. Do not destroy, bend stretch, twist or cut the power cord to avoid electric shock.
12. People who are sensitive to heat should not use the heat function of product.
13. *Recommended usage time for an individual is 15 minutes at a stretch.*
14. Please consult your doctor before using if you have the following situations
  - Receiving medical treatment or experiencing physical discomfort
  - Suffering from hypertension, heart disease, cerebral arteriosclerosis and cerebral thrombosis
  - Suffering from cancer
  - Pregnant or in a period
  - Suffering from osteoporosis
  - Using a pacemaker or implantable medical devices
  - Suffering from acute diseases
  - Minors
  - If you have recently had surgery on your legs
15. Do not fall asleep when using the product.
16. Do not use this product when you are drunk or feeling uncomfortable.
17. Do not use this product over thirty minutes in a row.
18. Do not use if the product's cloth cover has been damaged.

19. Do not use the product one hour before and after eating so as to avoid adverse effects.
20. If the product is used immediately after being transferred from a cold environment to a warm environment, the accumulation of water vapor inside the product may affect the normal operation of the product or even cause mechanical failure. You are suggested to put the product in a normal room temperature for 1 hour before starting to use it after the transfer.
21. Do not walk or do sport while using the leg massager

### **C. HOW TO USE IT**

#### **1) Battery charging**

Before first use , connect the USB cable to the charging port of the controller. Connect the power adaptor provided to a wall socket. The Power/charge indicator LED will turn red color during the charging state and will become green when the battery is fully charged.

Unplug the device once battery is fully charged.

You can use your massager 2H before next charging.

#### **2) Operating instructions**

- Check that the device is OFF before to install
- Do not apply or remove cuff when the chambers are inflated as it may cause damage to the product.
- Application of cuff is same for both legs.
- Wrap the cuff around the calf.
- Long press ON/OFF button for 3 seconds to start the device and the LED indicator will turn ON.
- The appliance switches directly to mode 1 (kneading), intensity 1 (low) and heat level 1 (low).
- Use the mode selection button ② to switch to "pressure" (= pressure) or "auto" (automatic mode). To select the desired program:
  - o Kneading = kneading mode
  - o Pressure = pressure mode
  - o Auto = automatic mode
- Use the intensity selection knob ③ to select one of the 3 intensities Low/Mid/High according to your preference.
- Use the heat setting knob ④ to select one of the 2 Low//High (= weak, strong) heat intensities according to your preferences.
- When the session ends, the device will go into standby mode. (around 20min / session)
- If you don't want to use the device for the time being and turn off the device, press the on/off ① button once.

Note :

When the chamber is inflating it is normal to hear a humming noise and feel a tightening sensation around your leg.

In the event of a pressure error or the chamber's pressure is unstable, all the LED on the device will turn on. Wait until the chambers fully deflated before pressing any button.

### **D. CARE AND MAINTENANCE**

Clean the product only with a dry or slightly damp, lint-free cloth.

Do not use abrasive cleaners to clean the unit.

### **E. FEATURES**

- Materials: 100% polyester fabric, ABS
- Power supply : Li-ion battery 2000mAh, 3.7V, 6W USB charging cable included (AC adapter not included) - Battery life : about 2h
- Massage: 3 modes & 3 intensity levels - Heat: 2 intensities

*The European directive 2012/19/EU about electronic and electric waste, requires that you can't throw away defective domestic appliance with common waste. Used device has to be collect separately to optimize the recuperation rate and the recycling of the materials to reduce the impact on the health and the environment.*

***Precision and precautions about the battery use***

M 7117 / FE-7117

- *Throwing away a used battery in the environment or with household waste pollutes and avoids the recuperation of reusable materials.*
- *It is important to limit the battery consumption and to respect the following rules:*
  - *Favor alkaline battery (last longer than saline battery) and when it's possible use rechargeable battery.*
  - *Leave used batteries and accumulators in adapted containers that you can find in some retailers. In that way, metals don't pollute the environment and the health.*
- *Battery has to be inserted by respecting the polarity show on the apparel and the battery. An incorrect positioning can damage the apparel, cause leaking and even provoke fire and the explosion of the battery.*
- *To assure a good functioning, battery has to be in a good condition. If there is a problem, put new battery.*
- *Never try to recharge battery not made for this. It could leak, warm up, cause fire or explode.*
- *Replace all the battery at the same time. Never mix up alkaline battery with saline battery.*
- *Used battery has to be removed.*
- *Also, remove the battery if you're not using the apparel for a long time.*
- *Never try to short-circuit the terminal of a battery.*
- *Never throw away the battery in fire because they can explode.*
- *Take out the accumulator from the apparel before reload it.*
- *We recommend that an adult watch if a child change the battery to respect those rules.*
- *If you swallow a battery, go immediately to your doctor or the hospital.*
- *For the best ongoing performance, we recommend using only fresh, high-quality, alkaline batteries.*

**Importé par / Importiert von / Geïmporteerd door / Imported by :**  
**CAPTELEC, 59170 Croix - FRANCE**

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